

Dispo Septembre 2016
 PVC: 39,99€



BF 195 Impédancemètre en verre BF 195 Diagnostic scale

Design attrayant et écran LCD arrondi

Ecran XL à rétro-éclairage bleu

Plateau moderne

Taille des chiffres 28 mm

Plateau de pesée en verre de sécurité

Electrodes en acier inoxydable

Affichage: poids, graisse corporelle, eau corporelle, taux de muscle, calories requises et IMC

Graduation:

- 0,1% pour la graisse et l'eau corporelles et le taux de muscle

- 100 g pour le poids et la masse osseuse

Sélecteur kg/lb/st

10 mémoires utilisateurs

5 niveaux d'activité

Démarrage et arrêt automatique

Capacité de charge 180 kg

Graduation 100 g

Indication de surcharge

1 pile x 3V CR 2032 incluse

Dimensions: 31 x 31 x 2,3 cm

Poids: env. 1670 g

Garantie 5 ans

UE: 4 / Carton de transport: -

EAN-Nr.: 4211125 74816 6

Art-Nr.: 748.16

Attractive, round LCD display

Blue illuminated XL display

Modern frame

Digit size 28 mm

Weighing area made from safety glass

Electrodes made of brushed stainless steel

Determination of body weight, fat, water, muscle density, bone mass, calorie requirement (AMR)

Introduction:

- 0.1% for body fat, body water and muscle density

- 100 g for body weight and bone mass

Change over to kg/lb/st

10 user memories

5 activity levels

Automatic switch-on and switch-off

Capacity 180 kg

Graduation 100 g

Overload indication

1 x 3V CR 2032 battery included

Product measurements: 31 x 31 x 2.3 cm

Product weight: approx. 1670 g

5 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 74816 6

Item no.: 748.16