



## HK 55 Coussin multi-fonctions HK 55 Easyfix heating pad

Forme ergonomique pour le ventre, le dos et les articulations

Utilisable comme coussin pour la nuque une fois roulé

Fixations rapides et simples grâce à une technologie affinée

Avec des extensions élastiques - pour toutes tailles

Fabriqué en velours respirant

Souple et agréable sur la peau

3 niveaux de température éclairés

Arrêt automatique après env. 90 min.

Chauffe rapidement

Contrôle électronique de température

Système de sécurité BSS

Interrupteur amovible

Lavable en machine à 30°C

Öko-Tex 100  
répondent aux exigences écologiques et rigoureuses de la certification Oeko-Tex Standard 100, et vérifié par l'Institut de recherches Hohenstein.

Dimensions: 59 x 30 cm

Poids: env. 410 g

100 watts

Garantie 3 ans

UE: 4 / Carton de transport:  
EAN-Nr.: 4211125 21413 5  
Art.-Nr.: 214.13



Ergonomic shape for abdomen, back and joints

Rolled up also usable as neck heating pad

Particularly easy fixing by refined attachment technology

With elastic extensions - suitable for any size

Made of breathable velours

Cosy and skinfriendly

3 illuminated temperature settings

Automatic switch off after approx. 90 minutes

Rapid heat up

Electronic temperature control

BSS (Beurer Safety System)

Detachable control unit

Machine washable at 30°C

Oeko-Tex 100  
The textiles used for this device meet the stringent human ecological requirements of Oeko-Tex Standard 100, as verified by Hohenstein Research Institute.

Product measurements: 59 x 30 cm

Product weight: approx. 410 g

100 watt

3 year guarantee

Sales unit: 4 / Shipping carton:  
EAN no.: 4211125 21413 5  
Item no.: 214.13

